- Low starting resistance (2.3 kg / 5lbs)
- Cable column adjusts to provide exercise variety with 35 positions
- 1:2 pulley ratio provides low incremental weight changes and moderate speed work
- MED CE compliant



\* Shown with optional incremental weights and medical handles.

ECH SPECS	
Product Dimensions (L x W x H)	108 x 125 x 228 cm / 43" x 50" x 90"
Product Weight	174 kg / 384 lbs.
Weight Stack	68 kg / 150 lbs.
Effective Resistance	34 kg / 75 lbs.
Incremental Weight (optional)	Qty 4, 0.9 kg / 2 lbs. (.45 kg / 1 lb. effective resistance)
May Training Weight	71.6 kg / 158 lbs

